

CAIRS-PRO


Advanced technology

CAIRS SUITE to empower and make your surgery easier.

Designed by Prof. C. Mazzotta M.D., PhD and M. Zagari M.D.

J2810.4A Double intrastromal dissector spatula for loading CAIRS Rings segment

J4377 CAIRS holder to size the ring segment



JANACH[®]
WHERE VISION MEETS PRECISION

CAIRS-PRO

Advanced technology

SUGGESTED SET:

CAIRS INSTRUMENTS SUITE:



J2810.4A Double intrastromal dissector spatula for loading CAIRS Rings segment

J4377 CAIRS holder, to size the ring segment

ADDITIONAL TOOLS:

J2031.13 Adjustable speculum, with gentle movement (peek insert), 13 mm open anatomic blades, titanium

J2115.18A Double-ended instrument to unstick and lift the corneal flap in LASIK procedures. The delicate hook works at the flap edges inside the groove. The "C" shaped tapered loop, once positioned below the flap, lifts it with a gentle push on the horizontal plane (from the hinge to the exit point)

JD2890.600 (Ø 6.00) - **JD2890.900** (Ø 9.00) Disposable cutting through vacuum donor punch

J2820.4 Caliper, straight, 20 mm spread, graduated scale with 0.5 mm stepst

MDP30 Depth Knife

SUTURING TOOLS:

J3277 Tying and stitch forceps for nylon and monofilament sutures, straight, with delicate 6 mm tying platform. Recommended for 8-0 to 10-0 sutures

J4197 SnagFree Twister needle holder, curved, without lock, delicate tips, overall length: 128 mm, titanium. Recommended for 10-0 to 11-0 sutures

JANACH[®]
WHERE VISION MEETS PRECISION

Via Negretti, 3 - 22100 Como - Italy
Tel. +39 031574088 - international@janach.it - janach.it

© Copyright: E.Janach
BRO.124.04.26.EN